



FRITZ RUCK
Ophthalmologische Systeme GmbH

CATALOGUE

2010

PentaSys accessories

PentaSys 2	2
Sytem accessories	3 - 6
General consumables	7 - 8
Diathermy	9 - 11
I/A handles	12 - 13
Phaco	14 - 17
Vitrectomy	18 - 24
Silicone oil infusion / aspiration	25 - 27
Light fibre	28
Tubes	29 - 31
Tube Sets	31 - 34

Consumables ophthalmic surgery

Viscoelastics	35
BSS	35
Spare tubes	36

Fritz Ruck Ophthalmologische Systeme GmbH • Jülicher Straße 115 • 52249 Eschweiler /Germany
Tel.: +49 (0)2403 / 9455 - 0 • Fax: +49 (0)2403 / 9455 - 5
eMail: info@fritzruck-gmbh.de • www.fritzruck-gmbh.de

We always use the biggest effort to the fact that all contained information is entire, correct and topical. Errors and omissions excepted. At last, hence, for completeness, correctness and the last actuality no guarantee can be taken over.

PentaSys 2



PentaSys2 DPS1

dual pump system with peristaltic and venturi pump



PentaSys2 VPS1

venturi pump system



PentaSys2 DPS2

dual pump system with peristaltic and venturi pump and fast venting system



PentaSys2 VPS2

venturi pump system with venturi pump, backflush pump and fast venting system

System accessories

top mat

article-no.: 01P201 PU: 1 pieces

P2



for PentaSys2

top mat

article-no.: 01PS01 PU: 1 pieces

P1



for PentaSys1

hood

article-no.: 01PS02 PU: 1 pieces

P1 P2



pen holder

article-no.: 01P230 PU: 1 pieces

P1 P2



holder for 2OM50

pen


article-no.: 02OM50 PU: 1 pieces

P1 P2




for sterile touch-screen usage

PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

System accessories

trolley

article-no.: 01P248 PU: 1 pieces



P2

base height 85 cm, 95 cm and 105 cm.

Article-no.:
1P248 (85 cm)
1P249 (95 cm)
1P250 (105 cm)

additional tray

article-no.: 01P270 PU: 1 pieces



P2

size: 55cm x 45cm, max load 15kg

holder for collection bag

article-no.: 01JS20 PU: 1 pieces



P1

stainless steel, additional hole must be drilled

arm for instrument tray

article-no.: 01P255 PU: 1 pieces



P2

instrument tray, big

article-no.: 01P260 PU: 1 pieces



P2

size: 51cm, 43cm

P1 = This article can be used with **PentaSys 1** PU = Packing Unit

P2 = This article can be used with **PentaSys 2** **STERILE** = sterile article

System accessories

instrument tray, small

article-no.: 01P261 PU: 1 pieces

P2



size: 41cm, 33cm

infusion pole, electric

article-no.: 01P237 PU: 1 pieces

P2



holder 1P237 required

infusion pole, manual

article-no.: 01P216 PU: 1 pieces

P2



holder 1P237 required

holder for infusion pole

article-no.: 01P238 PU: 1 pieces

P2



for manual and electric infusion pole

footswitch with 4 bumpers


article-no.: 01FS10 PU: 1 pieces

P1 P2




linear right

PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

System accessories

footswitch

article-no.: 01PS14

PU: 1 pieces



linear right

 = This article can be used with **PentaSys 1** **PU** = Packing Unit

 = This article can be used with **PentaSys 2** **STERILE** = sterile article

General consumables

cassette

article-no.: 03P220 PU: 10 pieces

 disposable article

 P2



for systems with fast-venting black cassette holder

cassette

article-no.: 03P221 PU: 10 pieces

 disposable article

 P1  P2



for systems with silver cassette holder

monitor cover

article-no.: 03OM51 PU: 50 pieces

 disposable article

 P1



for sterile touch-screen usage

collection bag

article-no.: 03OM58 PU: 50 pieces

 P1  P2



with tap and non-return valve, 2000 ml

external filter


article-no.: 03PS22 PU: 3 pieces

 P1  P2




liquid and bacillus protection

PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

General consumables

aspiration sensor

article-no.: 02PS20 PU: 1 pieces



P1

for PentaSys1 with reusable sensor

paper for data printer

article-no.: 03PS25 PU: 5 pieces



P1 P2

aspiration sensor, disp.

article-no.: 03PS50 PU: 20 pieces



 disposable article

P1 P2

for PentaSys1 and PentaSys2 DPS/PPS with sensor module / silver cassette holder

aspiration sensor for fast venting

article-no.: 03P250 PU: 1 pieces



 disposable article

P2

for PentaSys2 DPS with fast venting / black cassette holder

 = This article can be used with **PentaSys 1** PU = Packing Unit

 = This article can be used with **PentaSys 2**  = sterile article


Diathermy

diathermy pen "sclera"

article-no.: 02BI70

PU: 1 pieces

P1 **P2**




straight, angular, 2-pin connector

diathermy pen "sclera"

article-no.: 02BI71

PU: 1 pieces

P1 **P2**




sideways, 2-pin connector

diathermy pen "sclera"

article-no.: 02BI72

PU: 1 pieces

P1 **P2**




blunt, buckled, 2-pin connector

diathermy forceps

article-no.: 02BI60

PU: 1 pieces

P1 **P2**




straight, 2-pin connector

diathermy forceps

article-no.: 02BI61

PU: 1 pieces

P1 **P2**



buckled, 2-pin connector

PU = Packing Unit

P1 This article can be used with **PentaSys 1**

STERILE = sterile article

P2 This article can be used with **PentaSys 2**


Diathermy

diathermy pen "Endo", 20G

article-no.: 02B173

PU: 1 pieces

P1 P2




straight, blunt, 2-pin connector, for intraocular use

diathermy pen "Endo", 20G

article-no.: 02B174

PU: 1 pieces

P1 P2




straight, spiky, 2-pin connector, for intraocular use

diathermy pen "Endo", 23G

article-no.: 14-7007

PU: 1 pieces

P1 P2



straight, 2-pin connector

diathermy pen "Endo", 20G, disposable


article-no.: 14-5004

PU: 10 pieces

disposable article

STERILE

P1 P2



straight, 2-pin connector

diathermy pen "Endo", 23G, disposable


article-no.: 03B123

PU: 5 pieces

disposable article

STERILE

P1 P2



straight, 2-pin connector

= This article can be used with **PentaSys 1** **PU** = Packing Unit

= This article can be used with **PentaSys 2** **STERILE** = sterile article

Diathermy

diathermy pen "Endo", 25G, disposable


article-no.: 03BI25

PU: 5 pieces

 disposable article

STERILE



straight, 2-pin connector

diathermy cable

article-no.: 02BI44

PU: 1 pieces



with round metal connector

diathermy cable

article-no.: 02BI46

PU: 1 pieces



with 2-pin connector

diathermy adapter for light fibre 03LI43

article-no.: 02BI50


PU: 1 pieces




usage with 3LI43 light-fibre including diathermy function

PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

I/A handles

bimanual handle aspiration, 21 G

article-no.: 021A36

PU: 1 pieces

P1 P2




ø 0.8mm, port Ø 0.35 mm, roughened tip

bimanual handle aspiration, 21 G

article-no.: 021A37

PU: 1 pieces

P1 P2




ø 0.8mm, port Ø 0.25 mm, roughened tip

bimanual handle infusion, 21 G

article-no.: 021A38

PU: 1 pieces

P1 P2



ø 0.8mm, 2 ports Ø 0.45 mm, roughened tip

bimanual I/A-set 21G 0,25 port

article-no.: 031A52

PU: 50 pieces

 disposable article

P1 P2



consists of: ø 0,8mm aspiration port Ø 0.25 mm, infusion port Ø 0.45 mm

bimanual I/A-set 21G 0,35 port

article-no.: 031A53

PU: 50 pieces

 disposable article

P1 P2



consists of: ø 0,8mm aspiration port Ø 0.35 mm, infusion port Ø 0.45 mm

 = This article can be used with **PentaSys 1** **PU** = Packing Unit

 = This article can be used with **PentaSys 2** **STERILE** = sterile article


I/A handles

I/A-handle

article-no.: 02IA34

PU: 1 pieces

P1 P2




16G/19G, 30°, with sleeve

I/A-handle

article-no.: 02IA33

PU: 1 pieces

P1 P2




16G/19G, 45°, with sleeve

I/A-handle

article-no.: 02IA32

PU: 1 pieces

P1 P2




16G/19G, straight, with sleeve

I/A-handle, metal sleeve

article-no.: 02IA30

PU: 1 pieces

P1 P2



16G/19G, straight

I/A-handle, metal sleeve

article-no.: 02IA31

PU: 1 pieces

P1 P2



16G/19G,45°

PU = Packing Unit

P1 This article can be used with **PentaSys 1**

STERILE = sterile article

P2 This article can be used with **PentaSys 2**

Phaco

phaco handle

article-no.: 02PH42 PU: 1 pieces

P1 P2




with connector for PentaSys1, PentaSys2, ProPhaco

steri box

article-no.: 06ST13 PU: 1 pieces

P1 P2



height 38mm

tool for phaco tip, metal

article-no.: 02PH58 PU: 1 pieces

P1 P2



stainless steel

tool for phaco tip

article-no.: 03PH68 PU: 10 pieces

P1 P2



plastic

protection cap

article-no.: 03PH70 PU: 10 pieces

P1 P2



for mounted phaco tip during sterilization

P1 = This article can be used with **PentaSys 1** PU = Packing Unit

P2 = This article can be used with **PentaSys 2** **STERILE** = sterile article

cleaning adapter

article-no.: 021A40

PU: 10 pieces



adapter male Luer / male Luer, for I/A-handles and phaco handles

phaco tip 'turbo', 30°, low bubble 1.2

article-no.: 03PH55

PU: 1 pieces



phaco tip, 30°, 1.2/0.9 x 0.9 mm, incision: 2.65 bis 3.2 mm, sleeve 03PH80 (blue)

phaco tip 'turbo', 30°, low bubble 0.82

article-no.: 03PH59

PU: 1 pieces



small incision phaco tip, 30°, 0.82/0.74 x 0.5 mm, incision 2.2 - 2.4 mm / sleeve 3PH87 (clear) or 2.5 - 2.65 mm / sleeve 03PH86 (brown)

phaco tip 'turbo', 30°,Mini

article-no.: 03PH63

PU: 1 pieces



micro phaco tip 30°, 0,8/0,68 x 0,5 mm, incision 1.6 - 2.0 mm, sleeve CF 662350 (orange)

phaco tip classic 'E', 30°, low bubble, 1.2


article-no.: 03PH50

PU: 1 pieces




phaco tip straight, 30°, 1.2 x 0.9 mm, incision: 2.65 - 3.2 mm, suggested sleeve 03PH80 (blue)

PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

Phaco

phaco tip classic 'E', 45°, low bubble, 1.2

article-no.: 03PH51 PU: 1 pieces

P1 P2



phaco tip straight, 45°, 1.2 x 0.9 mm, incision: 2.65 - 3.2 mm, suggested sleeve 03PH80 (blue)

phaco tip classic 'E', 30°, low bubble 0.9

article-no.: 03PH61- PU: 1 pieces

P1 P2



for bimanual and coaxial phaco, 30°, 0.9 x 0.7 mm, incision 2.5 - 2,65 mm / bimanual 1.4 mm, sleeve 03PH86 (brown)

phaco tip 'E', 45°

article-no.: 03PH61- PU: 1 pieces

P1 P2



standard needle, 45°, 0,9 x 0,7 mm

phaco tip 'fragmentation/Endo', 30°,

article-no.: 03PH60 PU: 1 pieces

P1 P2



tip for endo phaco \varnothing 0.9 mm / \varnothing 0.7 mm

test chamber

article-no.: 03PH90 PU: 10 pieces

P1 P2



P1 = This article can be used with **PentaSys 1** PU = Packing Unit

P2 = This article can be used with **PentaSys 2** **STERILE** = sterile article

sleeve 19 Gauge, 'standard'

article-no.: 03PH80 PU: 10 pieces



blue, 2.75-3.2mm incision

sleeve 19 Gauge, 'Cool Flow'

article-no.: 03PH85 PU: 4 pieces



yellow, rigid, 2.75-3.2mm incision

sleeve 20 Gauge

article-no.: 03PH86 PU: 6 pieces



brown, 2.4-2.65mm incision

sleeve 21 Gauge

article-no.: 03PH87 PU: 6 pieces



clear, 2.2-2.4mm incision

sleeve 23 Gauge

article-no.: CF602350 PU: 6 pieces



orange, 1.6-2.0mm incision

PU = Packing Unit

 This article can be used with **PentaSys 1**
STERILE = sterile article

 This article can be used with **PentaSys 2**

Vitrectomy

pneumatic cutter 20 G 1200


article-no.: 03VI60

PU: 1 pieces

 disposable article

STERILE



PentaSys2 0-1200cuts/min, PentaSys1 0-800cuts/min

pneumatic cutter 20 G "high speed" 3000


article-no.: 03VI66

PU: 1 pieces

 disposable article

STERILE





PentaSys2 0-3000cuts/min, PentaSys1 0-1000cuts/min

pneumatic cutter 23 G 3000


article-no.: 03VI04

PU: 1 pieces

 disposable article

STERILE



PentaSys2 0-3000cuts/min, PentaSys1 0-1000cuts/min

pneumatic cutter "ECO" 20 G

article-no.: 03VI50

PU: 5 pieces

 disposable article

STERILE



800cuts/min

pneumatic cutter 20 G, reuseable

article-no.: 03VI70

PU: 1 pieces

STERILE



0-800cuts/min, reuse 5 times

 = This article can be used with **PentaSys 1** PU = Packing Unit

 = This article can be used with **PentaSys 2**  = sterile article

Vitrectomy

electric cutter 0-2500cuts/min

article-no.: 02SC02 PU: 1 pieces

P1 P2



PentaSys1+2, adjustable port size

electric cutter blade

article-no.: 02SC20 PU: 1 pieces

P1 P2



port size 0,45 mm, 20 G

tube for electric cutter

article-no.: 03SC33 PU: 10 pieces

P1 P2



infusion sleeve for pneumatic cutter

article-no.: 03V134 PU: 5 pieces

 disposable article

STERILE

P1 P2




for anterior vitrectomy, usage with 20G cutter

infusion port 23 G, anterior chamber infusion

article-no.: M2658 PU: 5 pieces


 disposable article

STERILE




4.0mm, for anterior Vitrectomy

PU = Packing Unit

 This article can be used with **PentaSys 1**

 = sterile article

 This article can be used with **PentaSys 2**

Vitrectomy

scissors 23 G

article-no.: VR-1513 PU: 1 pieces



135 degr.

scissors 23 G

article-no.: VR-1511 PU: 1 pieces



straight

delicated grasping forceps 23 G

article-no.: VR-1553 PU: 1 pieces



Eckardt forceps 23 G

article-no.: VR-1537 PU: 1 pieces



serrated forceps 23 G

article-no.: VR-1520 PU: 1 pieces



P1 = This article can be used with **PentaSys 1** **PU** = Packing Unit

P2 = This article can be used with **PentaSys 2** **STERILE** = sterile article

Vitrectomy

infusion port 20 G

article-no.: 02VI30

PU: 1 pieces

P1 P2



5.0mm

infusion port 20 G

article-no.: 02VI29

PU: 1 pieces

P1 P2



6.0mm

infusion port 20 G

article-no.: 02VI28

PU: 1 pieces

P1 P2



4.0mm

infusions port 20 G


article-no.: 03VI15

PU: 10 pieces

 disposable article

STERILE

P1 P2



5.0mm

backflush needle 20 G

article-no.: 5151

PU: 5 pieces


 disposable article

STERILE




LuerLock

PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

Vitrectomy

brush needle 20 G

article-no.: 5017

PU: 5 pieces



 disposable article

STERILE

LuerLock

backflush needle 23 G

article-no.: 052221

PU: 5 pieces



 disposable article

STERILE

LuerLock

brush needle 23 G

article-no.: 052241

PU: 5 pieces



 disposable article

STERILE

LuerLock

Dual Bore PFC canula 20 G

article-no.: 5191

PU: 5 pieces



 disposable article

STERILE

Dual Bore PFC canula 23 G

article-no.: 5203

PU: 5 pieces



 disposable article

STERILE

 P1 = This article can be used with PentaSys 1 PU = Packing Unit

 P2 = This article can be used with PentaSys 2  STERILE = sterile article

backflush handle


article-no.: 052210 PU: 1 pieces



LuerLock

silicone inlay

article-no.: 052250 PU: 10 pieces



for backflush handle

 disposable article

vitrectomy set 23G "spot" 1-step

article-no.: 03VI23 PU: 1 pieces




consists of:
light fibre "spot", cutter 23G 3000 cuts/min, trocar set 1-step

 disposable article

STERILE  

vitrectomy set 23G "wide angle" 1-step

article-no.: 03VI24 PU: 1 pieces




consists of:
light fibre (wide angled /shielded), pneumatic cutter 23G 3000 cuts/min,
trocar set 1-step


 disposable article

STERILE  

PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

Vitrectomy

trocar set 23 G 1-step

article-no.: 03VI11

PU: 1 pieces



STERILE



consists of: for one-step technique, including plugs, infusion line, template and forceps

trocar set 23 G 1-Step ready for silicone oil infusion

article-no.: 03VI11-2

PU: 1 pieces



STERILE



consists of: for one-step technique, including plugs, infusion line, template and forceps

trocar set 23 G 2-step with valves

article-no.: 03VI12

PU: 1 pieces



STERILE



consists of: for two-step technique, with self-closing valves, including infusion line, template and forceps

trocar set 23 G 2-Step with valves ready for silicone oil infusion

article-no.: 03VI12-2

PU: 1 pieces



STERILE



consists of: for two-step technique, with self-closing valves, including infusion line, template and forceps

= This article can be used with **PentaSys 1** **PU** = Packing Unit

= This article can be used with **PentaSys 2** **STERILE** = sterile article

Silicone oil infusion / aspiration

infusion port 20 G, silicone oil

article-no.: 02ÖL30 PU: 1 pieces

P1 P2



5.0mm

infusion port 20 G, silicone oil

article-no.: 02ÖL29 PU: 1 pieces

P1 P2



6.0mm

infusion port 20 G, silicone oil

article-no.: 02ÖL28 PU: 1 pieces

P1 P2



4.0mm


infusion port 20 G, silicone oil, disposable

article-no.: 03ÖL30 PU: 10 pieces

 disposable article

STERILE

P1 P2



5.0mm


pressure tube for silicone oil injection

article-no.: 02ÖL32 PU: 1 pieces


P1 P2



PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

Silicone oil infusion / aspiration

silicone oil infusion connector

article-no.: 02ÖL10 PU: 1 pieces



box for silicone oil injection, high pressure

article-no.: 02ÖL31 PU: 1 pieces



adapter for silicone oil aspiration

article-no.: 02ÖL40 PU: 1 pieces



usage with 03ÖL43

adapter for silicone oil injection

article-no.: 02ÖL33 PU: 1 pieces



female Luerlock

adapter for silicone oil injection

article-no.: 02ÖL34 PU: 1 pieces



male Luerlock

 = This article can be used with **PentaSys 1** **PU** = Packing Unit

 = This article can be used with **PentaSys 2** **STERILE** = sterile article

Silicone oil infusion / aspiration

syringe 10 ml

article-no.: 03ÖL42 PU: 100 pieces

 **disposable article**

STERILE



for infusion

syringe 20 ml

article-no.: 03ÖL43 PU: 60 pieces

 **disposable article**

STERILE



for aspiration

pressure tube for silicone oil injection, disposable

article-no.: 03ÖL32 PU: 10 pieces


 **disposable article**

STERILE


 



PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

Light fibre

light fibre 20 G, spot


article-no.: 03LI40

PU: 5 pieces

 disposable article

STERILE



light fibre 20 G, wide angled


article-no.: 03LI42

PU: 5 pieces

 disposable article

STERILE



light fibre 23 G, spot


article-no.: 03VI01

PU: 5 pieces

 disposable article

STERILE



light fibre 23 G, wide angled, shielded

article-no.: 5210

PU: 5 pieces

 disposable article

STERILE



light fibre 20 G with diathermy

article-no.: 03LI43

PU: 5 pieces

 disposable article

STERILE



usage with 02BI50

 = This article can be used with **PentaSys 1** PU = Packing Unit

 = This article can be used with **PentaSys 2** **STERILE** = sterile article

I/A-tube PentaSys, disposable

article-no.: 03PS38 PU: 20 pieces

 disposable article

STERILE





for PentaSys1/2 with peristaltic pump

I/A-tube PentaSys incl. infusion bottle kit

article-no.: 03PS70 PU: 20 pieces

 disposable article

STERILE





03PS38 incl infusion bottle kit

I/A-tube PentaSys, reusable

article-no.: 03PS39 PU: 20 pieces

STERILE





for PentaSys1/2 with peristaltic pump, reusable

I/A-tube PentaSys2 VPS2, disposable

article-no.: 03P238 PU: 20 pieces

 disposable article

STERILE





PVC for PentaSys2 VPS2

I/A-tube PentaSys VPS1, disposable

article-no.: 03OM48 PU: 20 pieces

 disposable article


STERILE






PVC, for PentaSys1 VPS, PentaSys2 VPS1

PU = Packing Unit

 This article can be used with **PentaSys 1**

 = sterile article

 This article can be used with **PentaSys 2**

Tubes

I/A-tube PentaSys VPS1, disposable

article-no.: 03OM49 PU: 20 pieces

 disposable article

STERILE

 P1



silicone, for PentaSys1 VPS, PentaSys2 VPS1

I/A-tube PentaSys VPS1, reuseable

article-no.: 03OM47 PU: 20 pieces

STERILE

 P1




silicone, for PentaSys1 VPS, PentaSys2 VPS1

peristaltic pump tube PentaSys1

article-no.: 03PS37 PU: 5 pieces

 disposable article

 P1



for PentaSys1 with reusable sensor

air set for vitrectomy

article-no.: 03LU65 PU: 20 pieces

 disposable article

STERILE

 P1  P2



incl. filter and 3-way-valve

3-way tap

article-no.: 03LU60 PU: 50 pieces

 disposable article

STERILE

 P1  P2



 P1 = This article can be used with **PentaSys 1** PU = Packing Unit

 P2 = This article can be used with **PentaSys 2**  STERILE = sterile article

Tubes

filter for air infusion

article-no.: 03LU50 PU: 25 pieces



 disposable article

STERILE  

IOP Set

article-no.: 03LU66 PU: 20 pieces



 disposable article

STERILE  

Tube Sets

I/A-tube set PentaSys DPS No.1

article-no.: 03PS40 PU: 20 pieces



 disposable article

STERILE  

consists of I/A-tube 03PS38, 1 sleeve 03PH80 and 1 test chamber 03PH90

I/A-tube set PentaSys DPS No.6

article-no.: 03PS402 PU: 20 pieces




 disposable article

STERILE  

consists of I/A-tube 03PS38, 2 sleeves 03PH80 and 1 test chamber 03PH90

PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

Tube Sets

I/A-tube set PentaSys DPS No.3

article-no.: 03PS42 PU: 20 pieces

 disposable article

STERILE  




consists of I/A-tube 03PS38 and sensor 03PS50

I/A-tube set PentaSys DPS No.7

article-no.: 03PS79 PU: 20 pieces

 disposable article

STERILE  



consists of I/A-tube 03PS38 and 22G bimanual I/A-handles 0.45/0.35 port

I/A-tube set PentaSys DPS No.2

article-no.: 03PS41 PU: 20 pieces

 disposable article

STERILE  



consists of I/A-tube 03PS38, 1 sleeve 03PH80 and 1 test chamber 03PH90 , phaco tip 03PH65

I/A-tube set PentaSys DPS No. 8

article-no.: 03PS45 PU: 1 pieces

 disposable article

STERILE 



2.2mm incision

consists of I/A-tube 03PS38 with phaco tip 06PH59, 1 sleeve 03PH87 and 1 test chamber 03PH90

I/A-tube set PentaSys DPS No.5

article-no.: 03PS47 PU: 20 pieces

 disposable article

STERILE 



consists of I/A-tube 03PS39 and peristaltic pump tube 03PS37

 = This article can be used with **PentaSys 1** PU = Packing Unit

 = This article can be used with **PentaSys 2**  = sterile article

Tube Sets

I/A-tube set PentaSys2 VPS2 No.50

article-no.: 03P240 PU: 20 pieces

 disposable article

STERILE





consists of I/A-tube 03P238, 1 sleeve 03PH80 and 1 test chamber 03PH90


I/A-tube set PentaSys2 VPS2 No.52

article-no.: 03P242 PU: 20 pieces

 disposable article

STERILE





consists of I/A-tube 03P238, 1 sleeve 03PH80, 1 test chamber 03PH90 and bimanual I/A-handles 03IA35/3IA38

I/A-tube set PentaSys2 VPS2 No.53

article-no.: 03P271 PU: 20 pieces

 disposable article

STERILE





consists of I/A-tube 03P238 with infusion bottle kit, 1 sleeve 03PH80 and 1 test chamber 03PH90

I/A-tube set PentaSys2 VPS2 No.51

article-no.: 03P241 PU: 20 pieces

 disposable article

STERILE





consists of I/A-tube 03P238, 1 sleeve 03PH80, 1 test chamber 03PH90 and 1 phaco tip 03PH65

I/A-tube set PentaSys2 VPS No.54

article-no.: 03P243 PU: 1 pieces

 disposable article

STERILE






2.2mm incision

consists of I/A-tube 03P238 with phaco tip 06PH59, 1 sleeve 03PH87 and 1 test chamber 03PH90

PU = Packing Unit

 This article can be used with **PentaSys 1**

 = sterile article

 This article can be used with **PentaSys 2**

Tube Sets

I/A-tube set PentaSys VPS 1 No.30

article-no.: 03OM50 PU: 20 pieces

 **disposable article**

STERILE

P1



consists of I/A-tube 03OM48, 1 sleeve 03PH80 and 1 test chamber 03PH90

I/A-tube set PentaSys VPS 1 No.31

article-no.: 03OM46 PU: 20 pieces

 **disposable article**

STERILE

P1



consists of I/A-tube 03OM49, 2 sleeve 03PH80 and 1 test chamber 03PH90

P1 = This article can be used with **PentaSys 1** **PU** = Packing Unit

P2 = This article can be used with **PentaSys 2** **STERILE** = sterile article

Viscoelastica

MediHeal® Sodium Hyaluron 1,2 %

article-no.: 07VH01 PU: 1 pieces



 disposable article

STERILE

MediHeal® Sodium Hyaluron 1,4 %

article-no.: 07VH02 PU: 1 pieces



 disposable article

STERILE

MediHeal® Sodium Hyaluron 1,8%

article-no.: 07VH00 PU: 1 pieces



 disposable article

STERILE

MediCell® 2% Hydroxypropylmethylzellulose

article-no.: 07VC00 PU: 1 pieces



 disposable article

STERILE

BSS

BSS, plastic bottle, 250ml

article-no.: 07BS25 PU: 10 pieces




BSS, plastic bottle, 500ml


article-no.: 07BS51 PU: 10 pieces



PU = Packing Unit

 This article can be used with **PentaSys 1**

STERILE = sterile article

 This article can be used with **PentaSys 2**

Spare tubes

silicone tube

article-no.: 06SS01 PU: 1 pieces



50 m, \varnothing 1,0 x \varnothing 3,0 mm

silicone tube

article-no.: 06SS02 PU: 1 pieces



50 m, \varnothing 2,0 x \varnothing 4,0 mm

silicone tube

article-no.: 06SS03 PU: 1 pieces



50 m, \varnothing 3,0 x \varnothing 5,0 mm

silicone tube

article-no.: 06SS04 PU: 1 pieces



50 m, \varnothing 1,5 x \varnothing 3,5 mm

P1 = This article can be used with **PentaSys 1** **PU** = Packing Unit

P2 = This article can be used with **PentaSys 2** **STERILE** = sterile article

